
LC3B_Tandem_Puncta_QuantificationV2

Table of Contents

.....	1
Ask for userInputs	2
Create image data store	2
Initialize data storage	2
Start processing images	2
Create subvolumes for each channel	3
3D controlled watershed	3
Find the vacuoles for each cell identified	11
Accumulate results	20
Return counts for vacuoles as table	20

Code requires bioformats MATLAB Toolbox and MATLAB version R2021a or later.

Creator

Created by: Joaquin Quintana (last modified: 05-18-2021)

Email: Joaquin.Quintana@Colorado.edu

LC3B_Tandem_Puncta_QuantificationV2 is a batch processing function used to identify LC3B vacuoles labeled with eGFP and mCherry signals and is designed for researchers investigating autophagy. Code expects the users data to be grayscale 16-bit volume data either in the form of multidimensional tiff stacks or ND2 files with three channels corresponding to the eGFP, mCherry and nuclear marker for mammalian cells. Data used to test this function was acquired using a confocal microscope as out-of-focus light will introduce issues with segmentation and other modes of image acquisition that do not minimize out-of-focus light will not produce quality results.

The file calls the functions:

userInputs imreadVolume Watershed2 FindVacuoles TabularCountVacuoles

Data for single cell analysis is reported as the number of eGFP and mCherry vacuoles counted for each cell and is printed to the command window along with the fileName the data was extracted from and the cellindex corresponding to which cell the vacuoles were found in with respect to the image used. This data is returned as a table to the command window and can also be viewed by calling or opening the structure Results_Table.

NOTE: Each cell line is different and each optical system is as well. Users should experimentally determine the parameters to use for their data. A helper function is provided with the code to help users access and identify parameters for their data called VolumeThresholdhelper.

NOTE: Users can use MATLABs publishing option to capture all the outputs (recommended) for all figures created by the subfunctions above or place a waitfor() call after each figure to observe the results of each figure created as the processing is occurring. See MATLABs documentation for publishing options and using waitfor(). See github account for an example of the published output.

```
%clear the command window, users workspace and close all figures
clc;
clear;
close all;
```

Ask for userInputs

```
[Input] = userInput();
```

Create image data store

```
imds = imageDatastore(Input.Indir,...
    'IncludeSubfolders',true,'FileExtensions',
    {'.tif','ND2'},'LabelSource',...
    'foldernames');
imds_NumerOfFiles = size(imds.Files(),1);
```

Initialize data storage

```
dat = struct ;
dat.FileName = [];
dat.SingleCell = [];
```

Start processing images

```
for ii = 1:imds_NumerOfFiles

    %pull images from data store one at a time by indexing
    img = b fopen(imds.Files{ii,:});

    %get file name
    [pathstr,name] = fileparts(imds.Files{ii,:});

    .....
    .....

    .....
    .....

    .....
    .....

    .....
    .....

    %let the user know which image is being processed. Print the name to
    command
    %window
    disp("Begin processing image: ")
    disp(name)
```

```
%unpack user inputs for passing to functions
Nuclear_channel = Input.Nuclear_channel;
eGFP_channel = Input.eGFP_channel;
mCherry_channel = Input.mCherry_channel;
eGFP_threshold = Input.threshG;
mCherry_threshold = Input.threshR;
MinvacVoxeGFP = Input.minGFPvacVolume;
MinvacVoxmCherry = Input.minRFPvacVolume;
MinNucVox = Input.minNucleusVolume;
Solidity = Input.Solidity;

Begin processing image:
High numbers of eGFP and mCherry

Begin processing image:
Low numbers of eGFP puncta and high mCherry Puncta (2)

Begin processing image:
Low numbers of eGFP puncta and high mCherry Puncta

Begin processing image:
Low numbers of eGFP puncta and low numbers of mCherry Puncta (Noisy
Test)
```


Create subvolumes for each channel


```
[r,g,b] =
imreadVolume(img,Nuclear_channel,eGFP_channel,mCherry_channel);
```


3D controlled watershed

```
[bw_stats] = Watershed2(r,g,b,name,MinNucVox,true);
```

Data for plots generated from - High numbers of eGFP and mCherry










Data for plots generated from - Low numbers of eGFP puncta and high mCherry Puncta (2)



LC3B_Tandem_Puncta_QuantificationV2






Data for plots generated from - Low numbers of eGFP puncta and high mCherry Puncta



LC3B_Tandem_Puncta_QuantificationV2






Data for plots generated from - Low numbers of eGFP puncta and low numbers of mCherry Puncta (Noisy Test)

LC3B_Tandem_Puncta_QuantificationV2









Find the vacuoles for each cell identified

```
[SingleCellStructure] =  
FindVaculoes(r,g,b,bw_stats,name,eGFP_threshold,...  
  
mCherry_threshold,MinvacVoxeGFP,MinvacVoxmCherry,MinNucVox,Solidity);
```



LC3B_Tandem_Puncta_QuantificationV2





LC3B_Tandem_Puncta_QuantificationV2





LC3B_Tandem_Puncta_QuantificationV2




LC3B_Tandem_Puncta_QuantificationV2





LC3B_Tandem_Puncta_QuantificationV2





LC3B_Tandem_Puncta_QuantificationV2




LC3B_Tandem_Puncta_QuantificationV2



LC3B_Tandem_Puncta_QuantificationV2





Accumulate results

```
dat(ii).FileName = name;
dat(ii).SingleCell = SingleCellStructure;

end
```

Return counts for vacuoles as table

```
close all;
[Results_Table] = TabularCountVacuoles(dat);

disp(Results_Table)
imshow('Results.png')
snapnow
```

		<i>fileName</i>
<i>CellIndex</i>	<i>eGFP_Vacuoles</i>	<i>mCherry_Vacuoles</i>

{ 'High numbers of eGFP and mCherry'	1	4	5
{ 'High numbers of eGFP and mCherry'	2	2	4
{ 'High numbers of eGFP and mCherry'	3	2	4
{ 'Low numbers of eGFP puncta and high mCherry Puncta (2)'	1	1	3
{ 'Low numbers of eGFP puncta and high mCherry Puncta (2)'	2	3	7

```
{ 'Low numbers of eGFP puncta and high mCherry Puncta (2)'
}      3      4      8
{ 'Low numbers of eGFP puncta and high mCherry Puncta (2)'
}      4      0      2
{ 'Low numbers of eGFP puncta and high mCherry Puncta'
}      1      0      1
{ 'Low numbers of eGFP puncta and high mCherry Puncta'
}      2      3      8
{ 'Low numbers of eGFP puncta and high mCherry Puncta'
}      3      4      9
{ 'Low numbers of eGFP puncta and high mCherry Puncta'
}      4      2      7
{ 'Low numbers of eGFP puncta and low numbers of mCherry Puncta
(Noisy Test)'}      1      0      0
{ 'Low numbers of eGFP puncta and low numbers of mCherry Puncta
(Noisy Test)'}      2      0      0
{ 'Low numbers of eGFP puncta and low numbers of mCherry Puncta
(Noisy Test)'}      3      0      0
{ 'Low numbers of eGFP puncta and low numbers of mCherry Puncta
(Noisy Test)'}      4      0      0
{ 'Low numbers of eGFP puncta and low numbers of mCherry Puncta
(Noisy Test)'}      5      0      2
{ 'Low numbers of eGFP puncta and low numbers of mCherry Puncta
(Noisy Test)'}      6      0      0
```

Figure 1

fileName	CellIndex	eGFP_Vacuoles	mCherry_Vacuoles
{ 'High numbers of eGFP and mCherry' }	1	4	5
{ 'High numbers of eGFP and mCherry' }	2	2	4
{ 'High numbers of eGFP and mCherry' }	3	2	4
{ 'Low numbers of eGFP puncta and high mCherry Puncta (2)' }	1	1	3
{ 'Low numbers of eGFP puncta and high mCherry Puncta (2)' }	2	3	7
{ 'Low numbers of eGFP puncta and high mCherry Puncta (2)' }	3	4	8
{ 'Low numbers of eGFP puncta and high mCherry Puncta (2)' }	4	0	2
{ 'Low numbers of eGFP puncta and high mCherry Puncta' }	1	0	1
{ 'Low numbers of eGFP puncta and high mCherry Puncta' }	2	3	8
{ 'Low numbers of eGFP puncta and high mCherry Puncta' }	3	4	9
{ 'Low numbers of eGFP puncta and high mCherry Puncta' }	4	2	7
{ 'Low numbers of eGFP puncta and low numbers of mCherry Puncta (Noisy Test)'} }	1	0	0
{ 'Low numbers of eGFP puncta and low numbers of mCherry Puncta (Noisy Test)'} }	2	0	0
{ 'Low numbers of eGFP puncta and low numbers of mCherry Puncta (Noisy Test)'} }	3	0	0
{ 'Low numbers of eGFP puncta and low numbers of mCherry Puncta (Noisy Test)'} }	4	0	0
{ 'Low numbers of eGFP puncta and low numbers of mCherry Puncta (Noisy Test)'} }	5	0	2
{ 'Low numbers of eGFP puncta and low numbers of mCherry Puncta (Noisy Test)'} }	6	0	0